

Recombinant Human Growth Arrest-Specific Protein 7,GAS-7,KIAA0394 (N-6His)



Catalog No: #AP74957

Orders: order@signalwayantibody.com

Package Size: #AP74957-1 10ug #AP74957-2 50ug #AP74957-3 500ug #AP74957-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Growth Arrest-Specific Protein 7,GAS-7,KIAA0394 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O60861
Uniprot	O60861
GeneID	8522;
Calculated MW	49.4kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMVPPPPGEESQTVILPPGWQSYLSPQGRRYYVNTTTNETTWERPSSSPGIP ASPGSHRSSLPPTVNGYHASGTPAHPPEAHMSVRKSTGDSQNLGSSSPSKKQSKENTITINCVTFFHPDPM PEQQLLKPTIEWSYCDYFWADKKDPQNGTVAGFELLQKQLKQKQKEMSEFIRERIKIEEDYAKNLAKLS QNSLASQEEGSLGEAWAQVKSLADEAEVHLKFSAKLHSEVEKPLMNFRENFKKDMKKCDHHIADLRKQLAS RYASVEKARKALTERQRDLEMKTQQLLEIKLSNKTEEDIKARRKSTQAGDDLRCVDLYNQAQSKWFEEMVT TTLELERLEVERVEMIRQHLCQYTLRHETDMFNQSTVEPVDQLLRKVDPKADRELWVREHKTGNIRPVDMEI
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 2mM DTT, 10% Glycerol, pH 8.8.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only