

Recombinant Human γ-Glutamylaminocyclotransferase,GGACT (N-6His)

Catalog No: #AP74960

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74960-1 10ug #AP74960-2 50ug #AP74960-3 500ug #AP74960-4 1mg

Description

Product Name	Recombinant Human γ-Glutamylaminocyclotransferase,GGACT (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9BVM4
Uniprot	Q9BVM4
GeneID	87769;
Calculated MW	19.5kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMALVFVYGT LKRGQPNHRVLRDGAHGSAAFRARGRTLEPYPLVIAGEHNIP WLLHLPGSGRLVEGEVYAVDERMLRFLDDFESCPALYQRTVLRVQLLEDRAPGAEPPAPTAVQCFVYSRAT FPPEWAQLPHHDSYDSEGP HGLRYNPRENR
Formulation	Supplied as a 0.2 μm filtered solution of 20mM TrisHCl, 100mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only