

Recombinant Human Geranylgeranyl Pyrophosphate Synthase, GGPS1 (N-6His)



Catalog No: #AP74961

Orders: order@signalwayantibody.com

Package Size: #AP74961-1 10ug #AP74961-2 50ug #AP74961-3 500ug #AP74961-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Geranylgeranyl Pyrophosphate Synthase, GGPS1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O95749
Uniprot	O95749
GeneID	9453;
Calculated MW	37kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMEKTQETVQRILLEPYKYLQLPGKQVVRTKLSQAFNHWLKVPEDKLQIIIIEVT EMLHNASLLIDDIEDNSKLRRGFPVAHSIYGIPSVINSANYVYFLGLEKVLTLDHPDAVKLFTRQLLELHQGGGL DIYWRDNYTCPTEEEYKAMVLQKTGGLFGLAVGLMQLFSDYKEDLKPLLNTLGLFFQIRDDYANLHSKEYSEN KSFCEDLTEGKFSFPTIHAIWSRPESTQVQNILRQRTEENIDIKKYCVHYLEDVGSFEYTRNTLKELEAKAYKQID ARGGNPELVALVKHLSKMFKEENE
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0 .
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only