

Recombinant Human Hemoglobin Subunit ζ ,HBAZ (N-6His)

Catalog No: #AP74963

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74963-1 10ug #AP74963-2 50ug #AP74963-3 500ug #AP74963-4 1mg

Description

Product Name	Recombinant Human Hemoglobin Subunit ζ ,HBAZ (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P02008
Uniprot	P02008
GeneID	3050;
Calculated MW	17.8kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMSLTKTERTIIVSMWAKISTQADTIGTETLERLFLSHPQTKTYFPDFLHPGSA QLRAHGSKVVAAVGDAVKSIDDIGGALSKLSELHAYILRVDPVNFKLLSHCLLVTLAARFPADFTAEHAHAWDK FLSVSSVLTEKYR
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM TrisHCl, 100mM NaCl, 2mM DTT, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only