

Recombinant Human Inositol Monophosphatase 1, IMPA1, IIMPase 2 (N-6His)

Catalog No: #AP74966

Orders: order@signalwayantibody.com

Package Size: #AP74966-1 10ug #AP74966-2 50ug #AP74966-3 500ug #AP74966-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Inositol Monophosphatase 1, IMPA1, IIMPase 2 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P29218
Uniprot	P29218
GeneID	3612;
Calculated MW	32.3kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMADPWQECMDYAVTLARQAGEVVCEAIKNEMNVMKSSPVDLVTATDQKV EKMLISSIKEKYPSSHSFIGEESVAAGEKSILTDNPTWIIDPIDGTTNFVHRFPFVAVSIGFAVNKKIEFGVVYSCVE GKMYTARKGKGAFNGQKLVQVSQQEDITKLLVTELGSSRTPETVRMVLNMEKLFIPVHGIRSVGTAAVN MCLVATGGADAYYEMGIHCWDVAGAGIIVTEAGGVLMDVTGGPFDLMSRRVIAANNRILAERIAKEIQVIPLQR DDED
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.25.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only