

# Recombinant Human L-Lactate Dehydrogenase A Chain,LDH-A,PIG19 (N-6His)



Catalog No: #AP74967

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP74967-1 10ug #AP74967-2 50ug #AP74967-3 500ug #AP74967-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human L-Lactate Dehydrogenase A Chain,LDH-A,PIG19 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P00338
Uniprot	P00338
GeneID	3939;
Calculated MW	38.8kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMATLKDQLIYNLLKEEQTPQNKITVVGAVGMACAISILMKDLADELALVDVI EDKCLKGEMMDLQHGSLFLRTPKIVSGKDYNVTANSKLVITAGARQQEGESRLNLVQRNVNIFKFIIPNVVKYSP NCKLLIVSNPVDILTIVAWKISGFPKNRVIGSGCNLDSARFRYLMGERLGVHPLSCHGWVLGEHGDSSVPVW SGMNVAGVSLKTLHPDLGTDKDKKEQWKEVHKQVVESAYEVIKLKGYSWAIGLSVADLAESIMKNLRRVHPVS TMIKGLYGIKDDVFLSVPCILGQNGISDLVKVTLTSEEEARLKKSAADTLWGIQKELQF
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 50% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only