

# Recombinant Human L-Lactate Dehydrogenase B Chain,LDH-B (N-6His)



Catalog No: #AP74968

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP74968-1 10ug #AP74968-2 50ug #AP74968-3 500ug #AP74968-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human L-Lactate Dehydrogenase B Chain,LDH-B (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P07195
Uniprot	P07195
GeneID	3945;
Calculated MW	38.8kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMATLKEKLIAPVAEEEEATVPNNKITVVGVGQVGMACAISILGKSLADELALVD VLEDKCLKGEMMDLQHGSFLQTPKIVADKDYSVTANSKIVVVTAGVRQQEGESRLNLVQRNVNVFKFIIPQIVK YSPDCIIIVVSNPVDILTIVTWKLSGLPKHRVIGSGCNLDSARFRYLMAEKLGHPSSCHGWILGEHGDSSVAV WSGVNVAGVSLQELNPEMGTDNDSENWKEVHKMVVESAYEVIKLKGTYNWAIGLSVADLIESMLKNLSRIHP VSTMVKGMYGIENEVFLSLPCILNARGLTSVINQKLDDEVAQLKKSADTLWDIQDKDL
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only