

Recombinant Human ZBTB17,Miz-1 (N-6His)

Catalog No: #AP74969



Package Size: #AP74969-1 10ug #AP74969-2 50ug #AP74969-3 500ug #AP74969-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZBTB17,Miz-1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q13105
Uniprot	Q13105
GeneID	7709;
Calculated MW	22.33kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMDFPQHSQHVLEQLNQQRQLGLLCDCTFVVDGVHFKAHKAVLAACSEYFK MLFVDQKDVVHLDISNAAGLGQVLEFMYTAKLSLSPENVDDVLAVATFLQMQDIITACHALKSLAEPATSPGGN AEALATEGGDKRAKEEKVATSTLSRLEQAGRSTPIGSRDLKEERGGQAQSAASGAEQTEKADA
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.25.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only