

Recombinant Human NHP2-Like Protein 1,NHP2L1 (N-6His)

Catalog No: #AP74972

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74972-1 10ug #AP74972-2 50ug #AP74972-3 500ug #AP74972-4 1mg

Description

Product Name	Recombinant Human NHP2-Like Protein 1,NHP2L1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P55769
Uniprot	P55769
GeneID	4809;
Calculated MW	16.3kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMTEADVNPAYPLADAHLLTKLLDLVQQSCNYQLRKGANEATKTLNRGISE FIVMAADAEPLEIILHLPLLCEDKNVPYVFRSKQALGRACGVSRPVIACSVTIKEGSQLKQIQSIQQSIERLLV
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 600mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only