

# Recombinant Human Nucleoside Diphosphate Kinase A,NDPKA (N-6His)



Catalog No: #AP74973

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP74973-1 10ug #AP74973-2 50ug #AP74973-3 500ug #AP74973-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Nucleoside Diphosphate Kinase A,NDPKA (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P15531
Uniprot	P15531
GeneID	4830;
Calculated MW	19.3kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMANCERTFIAIKPDGVQRGLVGEIHKRFEQKGFRLVGLKFMQASEDLLKEHYV DLKDRPFFAGLVKYMHS GPVVAMVWEGLNVVKTGRV MLGETNPADSKPGTIRGDFCIQVGRNIIHGSDSVES AEKEIGLWFHPEELVDYTSCAQNWIYE
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 1mM DTT, 10% Glycerol, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only