

# Recombinant Human Cyclin-Dependent Kinase 4 Inhibitor C,CDKN2C (N-6His)



Catalog No: #AP74975

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP74975-1 10ug #AP74975-2 50ug #AP74975-3 500ug #AP74975-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Cyclin-Dependent Kinase 4 Inhibitor C,CDKN2C (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P42773
Uniprot	P42773
GeneID	1031;
Calculated MW	20.3kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMAEPWGNELASAAARGDLEQLTSLQNNVNVNAQNGFGRTALQVMKLGNP EIARRLLLRGANPDLKDRTGFAVIHDAARAGFLDTLQTLLEFQADVNIEDNEGNLPLHLAAKEGHLRVVEFLVK HTASNVGHRNHKGDACDLARLYGRNEVVSLMQANGAGGATNLQ
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only