

# Recombinant Human PAF-AH subunit $\beta$ ,PAFAHB (C-6His)

Catalog No: #AP74976

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP74976-1 10ug #AP74976-2 50ug #AP74976-3 500ug #AP74976-4 1mg

## Description

Product Name	Recombinant Human PAF-AH subunit $\beta$ ,PAFAHB (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P68402
Uniprot	P68402
GeneID	5049;
Calculated MW	26.6kDa
Target Sequence	SQGDSNPAAIPHAIEDIQGDDRWMSQHNRFVLDCKDKEPDVLFVGDMSVQLMQQYEIWRELFSPHLALNFGI GGDTTRHVLWRLKNGELENIKPKVIVVWGTNNHENTAEVAGGIEAIVQLINTRQPQAKIIVLGLLPRGEKPNP LRQKNAKVNQLLKVSLPKLANVQLLDTDGGFVHSDGAISCHDMFDLHLTGGGYAKICKPLHELMQLEETPE EKQTTIAVEHHHHHH
Formulation	Supplied as a 0.2 $\mu$ m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only