

Recombinant Human Pterin-4- α -Carbinolamine Dehydratase,PHS,PCBD1 (N-6His)



Catalog No: #AP74978

Orders: order@signalwayantibody.com

Package Size: #AP74978-1 10ug #AP74978-2 50ug #AP74978-3 500ug #AP74978-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Pterin-4- α -Carbinolamine Dehydratase,PHS,PCBD1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P61457
Uniprot	P61457
GeneID	5092;
Calculated MW	14.2kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMAGKAHRLSAEERDQLLPNLRVAGWNELEGRDAIFKQFHFKDFNRAFGFMT RVALQAEKLDHHPEWVFNKVVHITLSTHECAGLSERDINLASFIEQVAVSMT
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM DTT, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only