

# Recombinant Human Phosphomevalonate Kinase,PMVK (N-6His)

Catalog No: #AP74980

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP74980-1 10ug #AP74980-2 50ug #AP74980-3 500ug #AP74980-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Phosphomevalonate Kinase,PMVK (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q15126
Uniprot	Q15126
GeneID	10654;
Calculated MW	24.2kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMAPLGGAPRLVLLFSGKRKSGKDFVTEALQSRLGADVCAVLRRLSGPLKEQY AQEHGLNFQRLLDTSTYKEAFRKMIRWGEEKRQADPGFFCRKIVEGISQPIWLVSDDRVRVSDIQWFREAYGA VTQTVRVVALEQSRQQRGWVFTPGVDDAESECGLDNFGDFDWVIENHGVEQRLEEQLLENLIEFIRSRL
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 1mM DTT, 10% Glycerol, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only