

Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase H,PPIH (N-6His)

Catalog No: #AP74985

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74985-1 10ug #AP74985-2 50ug #AP74985-3 500ug #AP74985-4 1mg

Description

Product Name	Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase H,PPIH (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O43447
Uniprot	O43447
GeneID	10465;
Calculated MW	21.4kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMAVANSSPVNPVFFDVSIGGQEVGRMKIELFADVVPKTAENFRQFCTGEF RKDGVPIGYKGSTFHRVIKDFMIQGGDFVNGDGTGVASIYRGPFADENFKLRHSAPGLLSMANS GPSTNGCQ FFITCSKCDWLDGKHVVFVKIIDGLLVMRKIENVPTGPNNKPKLPVVISQCGEM
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only