

Recombinant Human Triosephosphate Isomerase, TIM (N-6His)



Catalog No: #AP74992

Orders: order@signalwayantibody.com

Package Size: #AP74992-1 10ug #AP74992-2 50ug #AP74992-3 500ug #AP74992-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Triosephosphate Isomerase, TIM (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P60174
Uniprot	P60174
GeneID	7167;
Calculated MW	28.8kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMAPSRKFFVGGNWKMNQRKQSLGELIGTLNAAKVPADTEVVCAPPTAYIDF ARQKLDPKIAVAAQNCYKVTNGAFTGEISPGMIKDCGATWVVLGHSERRHVFGESEDELIGQKVAHALAEGLV IACIGEKLDEREAGITEKVVFEQTKVIADNVKDWKVVLAYEPVWAIGTGKTATPQQAQEVHEKLRGWLKSNV SDAVAQSTRIYGGSVTGATCKELASQPDVDGFLVGGASLKPEFVDIINAKQ
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 1mM DTT, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only