

Recombinant Human Tropomyosin α -1 Chain, TPM1 (N-6His)

Catalog No: #AP74993

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74993-1 10ug #AP74993-2 50ug #AP74993-3 500ug #AP74993-4 1mg

Description

Product Name	Recombinant Human Tropomyosin α -1 Chain, TPM1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09493
Uniprot	P09493
GeneID	7168;
Calculated MW	35kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMDAIKKKMQMLKLDKENALDRAEQAEADKKAEDRSKQLEDELVSLQKKLK GTEDELDKYSEALKDAQEKLELAEKKATDAEADVASLNRRRIQLVEEELDRAQERLATALQKLEEAKEKADESE RGMKVIESRAQKDEEKMEIQEIQLKEAKHIAEDADRKYEEVARKLVIIESDLERAEEERAELSEGQVRQLEEQLRI MDQTLKALMAAEDKYSQKEDRYEEEIKVLSDKLKEAETRAEFAERSVTKLEKSIDDLEDELYAQKLKYKAISEEL DHALNDMTSM
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM Tris, 1mM DTT, 100mM NaCl, 20% Glycerol, pH 8.0 .
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only