

Recombinant Human UDP-Glucose 4-Epimerase,GALE (N-6His)

Catalog No: #AP74995

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74995-1 10ug #AP74995-2 50ug #AP74995-3 500ug #AP74995-4 1mg

Description

Product Name	Recombinant Human UDP-Glucose 4-Epimerase,GALE (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q14376
Uniprot	Q14376
GeneID	2582;
Calculated MW	40.44kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMAEKVLVTGGAGYIGSHTVLELLEAGYLPVVIDNFHNAFRGGGSLPESLRRV QELTGRSVEFEEMDILDQGALQRLFKKYSFMAVIHFAGLKAVGESVQKPLDYRVNLTGTIQLLEIMKAHGAVKN LVFSSSATVYGNPQYLPLDEAHPTGGCTNPYGKSKFFIEEMIRDLCQADKTWNAVLLRYFNPTGAHASGCIGE DPQGIPNNLMPYVSQVAIGRREALNVFGNDYDTEGTVGRDYIHVVDLAKGHIAALRKLKEQCGCRIYNLGTG TGYSVLQMVQAMEKASGKKIPYKVVAREGDVAACYANPSLAQEELGWTAALGLDRMCEDLWRWQKQNP S GFGTQA
Formulation	Supplied as a 0.2 µm filtered solution of 50mM TrisHCl, 150mM NaCl, 2mM DTT, 1mM EDTA, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only