

Recombinant Human Chloride Intracellular Channel Protein 2, CLIC2 (N-6His)



Catalog No: #AP74998

Orders: order@signalwayantibody.com

Package Size: #AP74998-1 10ug #AP74998-2 50ug #AP74998-3 500ug #AP74998-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Chloride Intracellular Channel Protein 2, CLIC2 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O15247
Uniprot	O15247
GeneID	1193;
Calculated MW	30.5kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMSGLRPGTQVDPEIELFVKAGSDGESIGNCPFCQRLF MILWLKGVKFNVTTV DMTRKPEELKDLAPGTNPPFLVYNKELKTD FIKIEEFLEQTLAPPRYPHLSPKYKESFDVGCNLF AKFSAYIKNT QKEANKNFEKSLLEKFKRLDDYLNTPLLDEIDPDSAEPPVSRRLFLDGDQLTLADC SLLPKLNIIKVA AKKYRD FDIPAEFSGVWRYLHNAYAREEFTHTCPEDKEIENTYANVAKQKS
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 1mM DTT, 20% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only