

Recombinant Human Uroporphyrinogen-III Synthase,UROIIIIS (C-6His)

Catalog No: #AP75000

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75000-1 10ug #AP75000-2 50ug #AP75000-3 500ug #AP75000-4 1mg

Description

Product Name	Recombinant Human Uroporphyrinogen-III Synthase,UROIIIIS (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P10746
Uniprot	P10746
GeneID	7390;
Calculated MW	29.7kDa
Target Sequence	MKVLLLKDAKEDDCGQDPYIRELGLYGLEATLIPVLSFEFLSLPSFSEKLSHPEDYGGLIFTSBRAVEAAELCLE QNNKTEVWERSLKEKWNKSVYVGNATASLVSKIGLDTEGETCGNAEKLAEYICRESSALPLLFPNGNLKR EILPKALKDKGIAMESITVYQTVAHPGIQGNLNSYSSQQGVPASITFFSPSGLTYSLKHIQELSGDNIDQIKFAAIG PTTARALAAQGLPVSCTAESPTPQALATGIRKALQPHGCCLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only