

Recombinant Human GABA Receptor-Associated Protein-Like 1, GABARAPL1 (N-6His)



Catalog No: #AP75001

Orders: order@signalwayantibody.com

Package Size: #AP75001-1 10ug #AP75001-2 50ug #AP75001-3 500ug #AP75001-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GABA Receptor-Associated Protein-Like 1, GABARAPL1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9H0R8
Uniprot	Q9H0R8
GeneID	23710;
Calculated MW	16.2kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMKFQYKEDHPFEYRKKEGKIRKKYPDRVPVIVEKAPKARVPDLDKRKYLVPSDLTVGQFYFLIRKRIHLRPEDALFFFVNNTIPPTSATMGQLYEDNHEEDYFLYVAYSDESVDYVYGGK
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only