

Recombinant Human NAD Kinase,NADK (N-6His)

Catalog No: #AP75002



Package Size: #AP75002-1 10ug #AP75002-2 50ug #AP75002-3 500ug #AP75002-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NAD Kinase,NADK (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O95544
Uniprot	O95544
GeneID	65220;
Calculated MW	44.4kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMSLHGPCPVTTFGPKACVLQNPQTIMHIQDPASQRLTWNKSPKSVLVIKKMR DASLLQPFKELCTHLMENMIVYVEKKVLEDPAIASDESFGAVKKKFKCTFREDYDDISNQIDFIICLGGDGTLLYA SSLFQGSVPPVMAFHLGSLGFLTPFSFENFQSQVTQVIEGNAAVLRSRLKVRVVKELRGKKTAVHNGLGEK GSQAAGLDMDVGKQAMQYQVLNEVVIDRGPSSYLSNVDVYLDGHLITTVQGDGVIVSTPTGSTAYAAAAGAS MIHPNVPAIMITPICPHLSFRPIVVPAGVELKIMLSPEARNTAWVSFDGRKRQEIRHGDSISITTSCYPLPSICVR DPVSDWFESLAQCLHWNVRKKQAHFEFFFFFFFFEEG
Formulation	Supplied as a 0.2 µm filtered solution of 50mM TrisHCl, 150mM NaCl, 1mM DTT, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only