

Recombinant Human DCN1-Like Protein 1,DCUN1D1 (N-6His)

Catalog No: #AP75003

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75003-1 10ug #AP75003-2 50ug #AP75003-3 500ug #AP75003-4 1mg

Description

Product Name	Recombinant Human DCN1-Like Protein 1,DCUN1D1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96GG9
Uniprot	Q96GG9
GeneID	54165;
Calculated MW	32.3kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMNKLSQKDKVRQFMIFTQSSEKTAVSCLSQNDWKLDVATDNFFQNP ELY IRESVKGSLDRKKLEQLYNRYKDPQDENKIGIDGIQQFCDDLALDPASISVLIWKFRAATQCEFSKQEFMDG MTELGCDSEIKLKAQIPKMEQELKEPGRFKDFYQFTFNFAKNPGQKGLDLEMAIYWNLVLNGRFKFLDLWNK FLEHHKRSIPKDTWNLDFSTMIADDMSNYDEEGAWPVLIDDFVEFARPQIAGTKSTTV
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only