

# Recombinant Human Uroporphyrinogen Decarboxylase,UROD (N-6His)

Catalog No: #AP75005

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75005-1 10ug #AP75005-2 50ug #AP75005-3 500ug #AP75005-4 1mg

## Description

Product Name	Recombinant Human Uroporphyrinogen Decarboxylase,UROD (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P06132
Uniprot	P06132
GeneID	7389;
Calculated MW	43kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMEANGLGPQGFPELKNDTFLRAAWGEETDYTPVWCMRQAGRYLPEFRETR AAQDFSTCRSPEACCELTQPLRRFPLDAAIIFSDILVVPQALGMEVTMVPKGKPSFPEPLREEQDLERLRDP EVSASELGYVFQAITLTRQLRAGRVPLIGFAGAPWTLMTYMVEGGGSSTMAQAKRWLYQRPQASHQLLRILT DALVPYLVGQVVAGAQAALQLFESHAGHLGPQLFNKFALPYIRDVAKQVKARLREAGLAPVPMIIFAKDGHFALE ELAQAGYEVVGLDWTVAPKKARECVGKTVTLQVNLDPICALYASEEEEIGQLVKQMLDDFGPHRYIANLGHGLY PDMDPEHVGAFFVDAVHKHSRLLRQN
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 1mM DTT, 1mM EDTA, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only