

Recombinant Human L-Xylulose Reductase,DCXR (N-6His)

Catalog No: #AP75009

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75009-1 10ug #AP75009-2 50ug #AP75009-3 500ug #AP75009-4 1mg

Description

Product Name	Recombinant Human L-Xylulose Reductase,DCXR (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q7Z4W1
Uniprot	Q7Z4W1
GeneID	51181;
Calculated MW	28.1kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMELFLAGRRVLVTGAGKGIGRGTQALHATGARVVAVSRTQADLDSLUREC PGIEPVCVDLGDWEATERALGSVGPVDLLVNNAVALLQPFLEVTKEAFDRSFEVNLRAVIQVSQIVARGLIAR GVPGAIVNVSSQCSQRAVTNHSVYCSTKGALDMLTKVMALELGPBKIRVNAVNPVMTSMGQATWSDPHK AKTMLNRIPLGKFAEVEHVVNAILFLLSDRSGMTTGSTLPVEGGFWAC
Formulation	Supplied as a 0.2 µm filtered solution of 50mM Tris, 150mM NaCl, 1mM DTT, 30% Glycerol, 1mM DTT, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only