

Recombinant Human Hydroxyacid Oxidase 1,HAO1 (N-Trx, 6His)



Catalog No: #AP75011

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75011-1 10ug #AP75011-2 50ug #AP75011-3 500ug #AP75011-4 1mg

Description

Product Name	Recombinant Human Hydroxyacid Oxidase 1,HAO1 (N-Trx, 6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UJM8
Uniprot	Q9UJM8
GeneID	54363;
Calculated MW	58.6kDa
Target Sequence	MSDKIIHLTDDSFDTDLKADGAILVDFWAEWCGPCKMIAPILDEIADEYQGKLTVAKLNIDQNPGTAPKYGIRGI PTLLLFKNGEVAATKVGALSKGQLKEFLDANLAGSGSGHMHSHHHSSGLVPRGSGMKETAALKFERQHMD SPDLGTDDDDKAMADIGSMLPRLICINDYEQHAKSVLPKSIYDYRSGANDEETLADNIAAFSRWKLYPRMLRN VAETDLSTSVLGQRVSMPICVGATAMQRMHV DGE LATVRACQSLGTGMMLSSWATSSIEEVAEAGPEALR WLQLYIKDREVTKLVRQAEKMGYKAIFVTVDTPYLG NR LDDVRNRFKLPQLRMKNFETSTLSFSPEENFG DDSGLAAYVAKAIDPSISWEDIKWLRRLTSLPIVAKGILRGDDAREAVKHGLNGLVSNHGARQLDGVPATIDVL PEIVEAVEGKVEVFLDGGVRKGTDLKALALGAKAVFVGRPIVWGLAFQGEKGVQDVLEILKEEFRLAMALSG CQNVKVIDKTLVRKNPLAVSKI
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only