

Recombinant Human SGTA (N-6His)

Catalog No: #AP75013



Package Size: #AP75013-1 10ug #AP75013-2 50ug #AP75013-3 500ug #AP75013-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SGTA (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O43765
Uniprot	O43765
GeneID	6449;
Calculated MW	36.2kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMDNKKRLAYAIQFLHDQLRHGGLSSDAQESLEVAIQCLETAFGVTVEDSDL ALPQTLPEIFEAAATGKEMPQDLRSPARTPPSEEDSAEAERLKTGNEQMKVENFEAAVHFYGGKAIELNPANA VYFCNRAAAYSKLGNYAGAVQDCERAICIDPAYSKAYGRMGLALSSLNKHVEAVAYYKALELDPDNETYKSN LKIAELKLREAPSPTGGVGSFDIAGLLNPNPGFMSMASNLNPNPQIQQLMSGMISGGNNPLGTPGTSPSQNDLA SLIQAGQQFAQQMQQQNPELIEQLRSQIRSRTPSASNDDQQE
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only