

Recombinant Human Z-DNA Binding Protein 1,ZBP1 (C-6His)

Catalog No: #AP75014

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75014-1 10ug #AP75014-2 50ug #AP75014-3 500ug #AP75014-4 1mg

Description

Product Name	Recombinant Human Z-DNA Binding Protein 1,ZBP1 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9H171
Uniprot	Q9H171
GeneID	81030;
Calculated MW	17.5kDa
Target Sequence	MAQAPADPGREGHLEQRILQVLTEAGSPVKLAQLVKECQAPKRELNQVLYRMKKELKVSLTSPATWCLGGTD PEGEGPAELALSSPAKRPQQAATIPETPGPQFSQQREEDIYRFLKDNGPQRALVIAQALGMRTAKDVNRDLY RMKSVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only