

Recombinant Human Ethanolaminephosphotransferase 1,EPT1 (N-GST)

Catalog No: #AP75015

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75015-1 10ug #AP75015-2 50ug #AP75015-3 500ug #AP75015-4 1mg

Description

Product Name	Recombinant Human Ethanolaminephosphotransferase 1,EPT1 (N-GST)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9C0D9
Uniprot	Q9C0D9
GeneID	85465;
Calculated MW	32.3kDa
Target Sequence	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDVKLTQSMAIIRYI ADKHNMLGGCPKERAEISMLEGAVLDIRYGVSRIAYSKDFETLKVDFLSKLPEMLKMFEDRLCHKTYLNGDHV THPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIQIDKYLKSSKYIAWPLQGWQATFGGGDHPKSDL VPRGSPEFHMAGYEYVSPEQLAGFDKYKYSVDTNPLSLYVMHPFWNTIVKVFPWLAP
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only