

## Recombinant Human Galectin-1, LGALS1 (C-6His)

Catalog No: #AP75017



Package Size: #AP75017-1 10ug #AP75017-2 50ug #AP75017-3 500ug #AP75017-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Galectin-1, LGALS1 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09382
Uniprot	P09382
GeneID	3956;
Calculated MW	15.78kDa
Target Sequence	ACGLVASNLNLKPGECLRVRGEVAPDAKSFVLNLGKDSNNLCLHFNPRFNAHGDANTIVCNSKDGGAWGTE QREAVFPFQPGSVAEVCITFDQANLTVKLPDGYEFKFPNRLNLEAINYMAADGDFKIKCVAFDLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 10mM PB, 200mM NaCl, 2mM DTT, pH 7.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only