

Recombinant Human MOB-Like Protein Phocein,MOB4,MOBKL3 (N-6His)



Catalog No: #AP75018

Orders: order@signalwayantibody.com

Package Size: #AP75018-1 10ug #AP75018-2 50ug #AP75018-3 500ug #AP75018-4 1mg

Support: tech@signalwayantibody.com

Description

| | |
|-----------------|---|
| Product Name | Recombinant Human MOB-Like Protein Phocein,MOB4,MOBKL3 (N-6His) |
| Host Species | E.coli |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | Q9Y3A3 |
| Uniprot | Q9Y3A3 |
| GeneID | 25843; |
| Calculated MW | 28.2kDa |
| Target Sequence | MGSSHHHHHHSSGLVPRGSHMVAEGTAVLRRNRPGTKAQDFYNWPDESFDMDSTLAVQQYIQQNIRAD CSNIDKILEPPEGQDEGVWKYEHLRQFCLELNGLAVKLQSECHPDTCTQMTATEQWIFLCAAHKTPKECPAID YTRHTLDGAACLLNSNKYFPSRVSIKESSVAKLGSVCRRIRIFSHAYFHHRQIFDEYENETFLCHRFTKFMKY NLMSKDNLIVPILEEEVQNSVSGESEA |
| Formulation | Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM DTT, pH 8.0. |
| Storage | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months. |

Note: This product is for in vitro research use only