

Recombinant Human Ubiquitin-Conjugating Enzyme E2 S,UBE2S (N-GST)



Catalog No: #AP75021

Orders: order@signalwayantibody.com

Package Size: #AP75021-1 10ug #AP75021-2 50ug #AP75021-3 500ug #AP75021-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ubiquitin-Conjugating Enzyme E2 S,UBE2S (N-GST)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q16763
Uniprot	Q16763
GeneID	27338;
Calculated MW	50.13kDa
Target Sequence	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDVKLTQSMARIYI ADKHNMLGGCPKERAEISMLEGAVLDIRYGVSRIAYS KDFETLKVDFLSKLPEMLKMFEDRLCHKTYLNGDHV THPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIQIDKYLKSSKYIAWPLQGWQATFGGGDHPKSDL VPRGSMNSNVENLPPHIIRLVYKEVTTLTADPPDGIKVFPNEEDLTLQVTIEGPEGTPYAGGLFRMKLLLGKD FPASPPKGYFLTKIFHPNVGANGEICVNLKRDWTAELGIRHVLLTIKCLLIHPNPESALNEEAGRLLLENYEEYA ARARLLTEIHGGAGGPGSRAEAGRALASGTEASSTDPGAPGGP GGAEGPMAKKHAGERDKKLAAKKKTDDK RALRRL
Formulation	Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only