

Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase-Like 1, PPIase, PPIL1(N-6His)

Catalog No: #AP75023

Orders: order@signalwayantibody.com

Package Size: #AP75023-1 10ug #AP75023-2 50ug #AP75023-3 500ug #AP75023-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase-Like 1, PPIase, PPIL1(N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9Y3C6
Uniprot	Q9Y3C6
GeneID	51645;
Calculated MW	20.4kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMAAIPDSWQPPNVYLETSMGIIVLELYWKHAPKTCKNFAELARRGYNGTK FHRIIKDFMIQGGDPTGTGRGGASIYGKQFEDELHPDLKFTGAGILAMANAGPDTNGSQFFVTLAPTQWLDGK HTIFGRVCQGIGMVNRVGMVETNSQDRPVDDVKIIKAYPSG
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only