

Recombinant Human Aldo-Keto Reductase 1C2,AKR1C2

Catalog No: #AP75027

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75027-1 10ug #AP75027-2 50ug #AP75027-3 500ug #AP75027-4 1mg

Description

Product Name	Recombinant Human Aldo-Keto Reductase 1C2,AKR1C2
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P52895
Uniprot	P52895
GeneID	1646;
Calculated MW	36.74kDa
Target Sequence	MDSKYQCVKLNDGHFMPVLGFGTYAPAEVPSKALEAVKLAIEAGFHHIDSAHVYNNEEQVGLAIRSKIADGS VKREDIFYTSKLWSNSHRPELVRPALERSLKNLQLDYVDLYLIHFPVSVKPGEEVIPKDENGKILFDTVDLCATW EAMEKCKDAGLAKSIGVSNFNHRLLEMILNKPGLKYKPVCNQVECHPYFNQRKLLDFCKSKDIVLVAYSALGS HREEPWVDPNSPVLLEDVPLCALAKKHKRTPALIALRYQLQRGVVVLAKSYNEQRIRQNVQVFEFQLTSEEMK AIDGLNRNVRYLTLDIFAGPPNYPFSDY
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 1mM DTT, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only