

Recombinant Human Inhibin β C Chain,INHBC (N-6His)

Catalog No: #AP75028

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75028-1 10ug #AP75028-2 50ug #AP75028-3 500ug #AP75028-4 1mg

Description

Product Name	Recombinant Human Inhibin β C Chain,INHBC (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P55103
Uniprot	P55103
GeneID	3626;
Calculated MW	14.83kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMGIDCQGGSRMCCRQEFFVDFREIGWHDWIIQPEGYAMNFCIGQCPLHIAG MPGIAASFHTAVLNLLKANTAAGTTGGGSCCVPTARRPLSLLYDRDSNIVKTDIPDMVVEACGCS
Formulation	Lyophilized from a 0.2 μ m filtered solution of 4mM HCl, 1mM DTT.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only