

Recombinant Human Arylsulfatase A,ARSA (C-6His)

Catalog No: #AP75033



Package Size: #AP75033-1 10ug #AP75033-2 50ug #AP75033-3 500ug #AP75033-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Arylsulfatase A,ARSA (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P15289
Uniprot	P15289
GeneID	410;
Calculated MW	52.9kDa
Target Sequence	RPPNIVLIFADDLGYGDLGCGYGHPSSTTPNLDQLAAGGLRFTDFYVPVSLCTPSRAALLTGRLPVRMGMYPGV LVPSSRGGLPLEEVTVAEVLAAARGYLTMAGKWHLGVGPEGAFPPHQGFHRFLGIPYSHDQGPCQNLTCFP PATPCDGGCDQGLVPIPLLANSVEAQPPWLPGLEARYMAFAHDLMAQAQRQDRPFFLYASHHHTYPQFS GQSFAERSGRGPFGDLSMELDAAVGTLMTAIGDLGLLEETLVIFTADNGPETMRMSRGGCSGLLRGCGKTTY EGGVREPALAFWPGHIAPGVTHELASSDLLPTLAALAGAPLPNVTLDGFDLSPLLLGTGKSPRQSLFFYPSTYP DEVRGVFAVRSGKYKAHFFTQGSASDATTADPACHASSSLTAHEPPLLYDLSKDPGENYNLLGGVAGATPEV LQALKQLQLLKAQLDAAVTFGPSQVARGEDPALQICCHPGCTPRPACCHCPDPAVDHDDHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only