

Recombinant Human Signaling Lymphocytic Activation Molecule,SLAM,CD150 (C-6His)

Catalog No: #AP75041

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75041-1 10ug #AP75041-2 50ug #AP75041-3 500ug #AP75041-4 1mg

Description

Product Name	Recombinant Human Signaling Lymphocytic Activation Molecule,SLAM,CD150 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q13291
Uniprot	Q13291
GeneID	6504;
Calculated MW	25.33kDa
Target Sequence	ASYGTGGRMMNCPKILRQLGSKVLLPLTYERINKSMNKSIHIVVTMAKSLNSVENKIVSLDPSEAGPPRYLGD RYKFYLENLTLGIRESRKEDEGWYLMTLEKNVSVQRFLQLRLYEQVSTPEIKVLNKTQENGTCCTLILGCTVEK GDHVAYSWSEKAGTHPLNPANSSHLLSLTLGPQHADNIYICTVSNPISNNSQTFSPWPGCRTDPSETKPDVHH HHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only