

## Recombinant Human OX-2,MOX1,CD200 (C-6His)

Catalog No: #AP75043



Package Size: #AP75043-1 10ug #AP75043-2 50ug #AP75043-3 500ug #AP75043-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human OX-2,MOX1,CD200 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P41217
Uniprot	P41217
GeneID	4345;
Calculated MW	23.48kDa
Target Sequence	QVQVVTQDEREQLYTPASLKCSLQNAQEALIVTWQKKKAVSPENMVTFSENHGVVVQIPAYKDKINITQLGLQN STITFWNITLEDEGCYMCLFNTFGFGKISGTACLTVYVQPIVSLHYKFSEDHLNITCSATARPAPMVFWKVP GIENSTVTLSHPNGTTSVTSILHIKDPKNQVGKEVICQLHLGTVTDFKQTVNKGVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only