

Recombinant Human Fc γ RIIa, FCGR2A, CD32a (C-6His)

Catalog No: #AP75047

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75047-1 10ug #AP75047-2 50ug #AP75047-3 500ug #AP75047-4 1mg

Description

Product Name	Recombinant Human Fc γ RIIa, FCGR2A, CD32a (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P12318
Uniprot	P12318
GeneID	2212;
Calculated MW	21.61kDa
Target Sequence	APPKAVLKLEPPWINVLQEDSVLTCQGARSPESDSIQWFHNGNLIPTHTQPSYRFKANNNDSGEYTCQTGQ TSLSDPVHLTVLSEWLVLQTPHLEFQEGETIMLRCHSWKDKPLVKVTFQNGKSQKFSRLDPTFSIPQANHSH SGDYHCTGNIGYTLFSSKPVTITVQVPSMGSSSPMGIIVDHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only