

Recombinant Human CD47,IAP,OA3 (C-6His)

Catalog No: #AP75050



Package Size: #AP75050-1 10ug #AP75050-2 50ug #AP75050-3 500ug #AP75050-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CD47,IAP,OA3 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q08722
Uniprot	Q08722
GeneID	961;
Calculated MW	14.76kDa
Target Sequence	QLLFNKTKSVEFTFCNDTVVIPCFVTNMEAQNTTEVYVKWKFKGRDIYTFDGALNKSTVPTDFSSAKIEVSQLL KGDASLKMDKSDAVSHTGNYTCEVTELTREGETIIEELKYRVVSWFSPVDHHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 10mM Tris-Citrate, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only