

# Recombinant Human Complement Factor B,CFB (C-6His)

Catalog No: #AP75056

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75056-1 10ug #AP75056-2 50ug #AP75056-3 500ug #AP75056-4 1mg

## Description

Product Name	Recombinant Human Complement Factor B,CFB (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P00751
Uniprot	P00751
GeneID	629;
Calculated MW	84.07kDa
Target Sequence	TPWSLAWPQGSCSLEGVEIKGGSFRLLEQEQALEYVCPSPGFYPYPVQTRTCRSTGWSWTLKTQDQKTVRKA ECRAIHCPRPDFENGGEYWPRSPYYNVSDIEISFHCYDGYTLRGSANRTCQVNGRWGQTAICDNGAGYCSN PGIPIGTRKVGSGYRLEDSVTYHCSRGLTLRGSQRRTCQEGGSSWGTPEPSCQDSFMYDTPQEVAEAFSSLT ETIEGVDAEDGHGPGEQQKRKIVLDPSGSMNIYLVLDGSDSIGASNFTGAKKCLVNLIEKVASYGVKPRYGLVT YATYPKIWVKVSEADSSNADWVTKQLNEINYEDHKLKSGTNTKKALQAVYSMMSWPDDVPPEGWNRTRHVII LMTDGLHNMGGDPITVIDEIRDLLYIGKDRKNPREDYLDVYVFGVGPLVNQVNINALASKKDNEQHVFKVKDM ENLEDVIFYQMIDESQSLSLCGMVWEHRKGTDYHKQPWQAKISVIRPSKGHESCMGAVVSEYFVLTAHCFTV DDKEHSIKVSVGGEKRDLEIEVVLPHPNYNINGKKEAGIPEFYDYDVALIKLKNLKYGQTIRPICLPCTEGTTRA LRLPPTTTCQQQKEELLPAQDIKALFVSEEEKLTRKEYVIKNGDKKGCERDAQYAPGYDKVKDISEVVTPRF LCTGGVSPYADPNTCRGDSGGPLIVHKRSRFIQVGVISWGVVDVCKNQKRQKQVPAHARDFHINLFQVLPWL KEKLQDEDLGFLVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only