

Recombinant Human CRTAMV,CD355 (C-6His)

Catalog No: #AP75057



Package Size: #AP75057-1 10ug #AP75057-2 50ug #AP75057-3 500ug #AP75057-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CRTAMV,CD355 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O95727
Uniprot	O95727
GeneID	56253;
Calculated MW	30.99kDa
Target Sequence	SLTNHTETITVEEGQTLTKCVTSLRKNSLQWLTPSGFTIFLNEYVPLKNSKYQLLHHSANQLSITVPNVTLQD EGVYKCLHYSDSVSTKEVKVIVLATPFKPILEASVIRKQNGEEHVLMCSTMRSKPPPQITWLLGNSMEVSGGT LHEFETDGKKCNTTSTLIIHTYGKNSTVDCIIRHRGLQGRKLVAPFRFEDLVTEETASDALERNLSLSSQDPQQ PTSTVSVTEDSSTSEIDKEEKEQTTQDPDLTTEANPQYLGLARKKSVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only