

Recombinant Human Protein FAM3D (C-6His)

Catalog No: #AP75068



Package Size: #AP75068-1 10ug #AP75068-2 50ug #AP75068-3 500ug #AP75068-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Protein FAM3D (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96BQ1
Uniprot	Q96BQ1
GeneID	131177;
Calculated MW	23.12kDa
Target Sequence	YMSFSMKTIRLPRWLAASPTKEIQVKKYKGLIKPCPANYFAFKICSGAANVVGPTMCFEDRMIMSPVKNNVG RGLNIALVNGTTGAVLGQKAFDMYSGDVMHLVKFLKEIPGGALVLVASYDDPGTKMNDESRLFSDLGSSYAK QLGFRDSWVFIGAKDLRGKSPFEQFLKNSPDTNKYEGWPELLEMEGCMPPKPFVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only