

Recombinant Human Fibronectin Leucine Rich Transmembrane Protein 1,FLRT1 (C-6His)



Catalog No: #AP75071

Orders: order@signalwayantibody.com

Package Size: #AP75071-1 10ug #AP75071-2 50ug #AP75071-3 500ug #AP75071-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Fibronectin Leucine Rich Transmembrane Protein 1,FLRT1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NZU1
Uniprot	Q9NZU1
GeneID	23769;
Calculated MW	56.52kDa
Target Sequence	IDSTTCPSVCRCDNNGFIYCNDRGLTSIPADIPDDATTLYLQNNQINNAGIPQDLKTKVNVQVIYLYENDLDEFPIN LPRSLRELHLQDNNVRTIARDSLARIPLEKLHLDDNSVSTVSIEEDAFADSKQLKLLFLSRNHLSSIPSGLPHTL EELRLDDNRISTIPLHAFKGLNSLRRLVLDGNLLANQRIADDTFSRLQNLTELSLVRNSLAAPPLNLP SAHLQKL YLQDNAISHIPYNTLAKMRELERLDLSNNNLTLPRGLFDDLGNLAQLLLRNNPWFQCGCNLMWLRDWWKARA AVVNVVRGLMCQGPEKVRGMAIKDITSEMDECFFETGPQGGVANAAKTTASNHASATTPQGSFLTKAKRPGL RLPDSNIDYPMATGDGAKTLAIHVKALTADSIRITWKATLPASSFRLSWLRRLGHSPAVGSITETLVQGDKTEYLL TALEPKSTYIICMVTMETSNAAYVADETPVCAKAETADSYGPTTTLNQEQNAGPMASLPVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only