

# Recombinant Human N-Acetylglucosamine-6-Sulfatase, GNS (C-6His)

Catalog No: #AP75074

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75074-1 10ug #AP75074-2 50ug #AP75074-3 500ug #AP75074-4 1mg

## Description

Product Name	Recombinant Human N-Acetylglucosamine-6-Sulfatase, GNS (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P15586
Uniprot	P15586
GeneID	2799;
Calculated MW	59.35kDa
Target Sequence	VFGVAAGTRRPNVVLLLTDDQDEVLGGMTPLKKTALIGEMGMTFSSAYVPSALCCPSRASILTGKYPHNHHV VNNTLEGNCSSKSWQKIQEPNTFPAILRSMCGYQTF FAGKYLNEYGAPDAGGLEHVPLGWSYWYALEKNSK YYNYTLSINGKARKKHGENYSVDYLTDLANVSLDFLDYKSNFEPFFMMIATPAPHSPWTAAPQYQKAFQNVFA PRNKNFNIHG TNKHWLIRQAKTPMTNSSIQFLDNAFRKRWQTLLSVDDLVEKLVKLEFTGELNNTYIFYTSDN GYHTGQFSLPIDKRQLYEFDIKVP LLVRGPGIKPNQTSKMLVANIDLGP TILDIAGYDLNKTQMDGMSLLPILRG ASNLTWRS DVLVEYQGEGRNVTDP TCPSLSPGV SQCFPDCV CEDAYNNTYACVRTMSALWNLQYCEFDDQ EVFVEVYNLTADPDQITNIAKTIDPELLGKMNYRLMMLQSCSGPTCRTPGVFDPGYRFDPRMLMFSNRGSRTR RFSKHL LDH HHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only