

Recombinant Human GZMB (C-6His)

Catalog No: #AP75075



Package Size: #AP75075-1 10ug #AP75075-2 50ug #AP75075-3 500ug #AP75075-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GZMB (C-6His)
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Accession No.	P10144
Uniprot	P10144
GeneID	3002
Calculated MW	26.6 KDa
Target Sequence	IIGGHEAKPHSRPYMAYLMIWDQKSLKRCGGFLIRDDFVLTAHCWGSSINVTLGAHNIKEQEPTQQFIPVKRPI PHPAYNPKNFSNDIMLLQLERKAKRTRAVQPLRLPSNKAQVKPGQTCSVAGWGQTAPLGKHSHTLQEVKMT VQEDRKCESDLRHYYDSTIELCVGDPEIKKTSFKGDSGGPLVCNKVAQGIVSYGRNNGMPPRACTKVSSFVH WIKKTMKRY
Formulation	Supplied as a 0.2 μ M filtered solution of PBS, pH7.4.
Storage	Store at β -70°C, stable for 6 months after receipt. Store at β -70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only