

Recombinant Human T Cell Immunoglobulin and Mucin Domain-3,TIM3,HAVCR2 (C-6His)



Catalog No: #AP75076

Orders: order@signalwayantibody.com

Package Size: #AP75076-1 10ug #AP75076-2 50ug #AP75076-3 500ug #AP75076-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human T Cell Immunoglobulin and Mucin Domain-3,TIM3,HAVCR2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8TDQ0
Uniprot	Q8TDQ0
GeneID	84868;
Calculated MW	20.86kDa
Target Sequence	SEVEYRAEVGQNAYLPCFYTPAAPGNLVPVCWGKGACPVFECGNVLRDTERDVNYWTSRYWLNDFRKG DVSLTIENVTLADSGIYCCRIQIPGIMNDEKFNLKLVIKPAKVTPAPTLQRDFTAAFPRLTTRGHGPAETQTLG SLPDINLTQISTLANELRDSRLANDLRDSGATIRVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only