

Recombinant Human High Mobility Group Protein B1, HMGB1 (C-6His)

Catalog No: #AP75077

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75077-1 10ug #AP75077-2 50ug #AP75077-3 500ug #AP75077-4 1mg

Description

Product Name	Recombinant Human High Mobility Group Protein B1, HMGB1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09429
Uniprot	P09429
GeneID	3146;
Calculated MW	25.93kDa
Target Sequence	GKGDPPKKPRGKMSSYAFFVQTCREEHKKKHPDASVNFSEFSKKCSERWKTMSAKEKGFEDMAKADKARY EREMKTYIPPKGETKKKFKDPNAPKRPPSAFFLFCSEYRPKIKGEHPGLSIGDVAKKLGEMWNNTAADDKQPY EKKAALKKEKYEKDIAAYRAKGPDAAKKGVVKAESKSKKKEEEDEEDEDEDEDEDEDEDEDEDEDEDEDEDEDEDEDDDEV DHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, 1mM DTT, 10% Glycerol, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only