

Recombinant Human Kallikrein 11, KLK11 (C-6His)

Catalog No: #AP75081



Package Size: #AP75081-1 10ug #AP75081-2 50ug #AP75081-3 500ug #AP75081-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Kallikrein 11, KLK11 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UBX7
Uniprot	Q9UBX7
GeneID	11012;
Calculated MW	26.28kDa
Target Sequence	ETRIIKGFECKPHSQPWQAALFEKTRLLCGATLIAPRWLLTAAHCLKPRYIVHLGQHNLQKEEGCEQTRTATES FPHPGFNNSLPNKDHRNDIMLVKMASPVSIWAVRPLTLSSRCVTAGTSLISGWGSTSSPQLRPLPHTLRCANI TIIHQKCEAYPGNITDTMVCASVQEGGKDCSCQGDSGGPLVCNQLQGIISWGDPCAITRKPVYTKVCKY VDWIQETMKNNVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 2mM CaCl ₂ , pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only