

## Recombinant Human Kallikrein 13, KLK13 (C-6His)

Catalog No: #AP75082



Package Size: #AP75082-1 10ug #AP75082-2 50ug #AP75082-3 500ug #AP75082-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Kallikrein 13, KLK13 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UKR3
Uniprot	Q9UKR3
GeneID	26085;
Calculated MW	28.03kDa
Target Sequence	GGVSQESSKVLNTNGTSGFLPGGYTCFPHSQPWQAALLVQGRLLCGGVLVHPKWVLTAAHCLKEGLKVYLG KHALGRVEAGEQVREVVHSIPHPEYRRSPHTLNHDHDIMLLELQSPVQLTGYIQTLP LSHNNRLTPGTTTCRVS GWGTTTSPQVNYPKTLQCANIQLRSDEECRQVYPGKITDNMLCAGTKEGGKDSCEGDSGGPLVCNRTLYGIV SWGDFPCGQPDRPGVYTRVSRVYLWIRETIVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, 10% Glycerol, pH 5.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only